

Product Datasheet

ATP1A3 antibody (HRP) (orb51239)

Description	ATP1A3 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Sodium/potassium-transporting ATPase subunit alpha-3 protein (143-278AA)
Target	ATP1A3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P13637
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.