

Product Datasheet

Allergen Bla g 4 antibody (orb51250)

Description	Allergen Bla g 4 antibody
Species/Host	Rabbit
Reactivity	Insect
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Blattella germanica Allergen Bla g 4 protein (13-182AA)
Target	Allergen Bla g 4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P54962
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.