

Product Datasheet

Major pollen allergen Aln g 1 antibody (FITC) (orb51256)

Description	Major pollen allergen Aln g 1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	FITC
Immunogen	Recombinant Alnus glutinosa Major pollen allergen Aln g 1 protein (2-160AA)
Target	Major pollen allergen Aln g 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P38948
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.