



Product Datasheet

EGFR (phospho-Tyr1092) antibody (orb5126)

Description	EGFR (phospho-Tyr1092) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human EGFR around the phosphorylation site of Tyr1092 (PE(p-Y)IN)
Target	phospho-EGFR (Tyr1092)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	170 kDa
Uniprot ID	P00533
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lung tissue.



Immunohistochemical staining of human lung tissue.