

Product Datasheet

cyp1a antibody (FITC) (orb51268)

Description	cyp1a antibody (FITC)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	FITC
Immunogen	Recombinant Zebrafish Cytochrome P450 protein (1-519AA)
Target	cyp1a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8UW07
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.