

Product Datasheet

nkx2.5 antibody (Biotin) (orb51273)

Description	nkx2.5 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Nkx2.5 protein (1-312AA)
Target	nkx2.5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q567L4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.