

Product Datasheet

CCDC74A antibody (orb513132)

Description	Rabbit polyclonal antibody to CCDC74A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 101-320aa of human CCDC74A protein
Target	CCDC74A
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	35kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.