

Product Datasheet

CFAP77 antibody (orb513209)

Description	Rabbit polyclonal antibody to CFAP77
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 1-173aa of human CFAP77 protein
Target	CFAP77
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	44kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.