

Product Datasheet

mrkA antibody (Biotin) (orb51321)

Description	mrkA antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Klebsiella pneumoniae Fimbrial subunit type 3 protein (23-202AA)
Target	mrkA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P12267
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.