

## Product Datasheet

### groL antibody (orb51342)

|                            |                                                                      |
|----------------------------|----------------------------------------------------------------------|
| <b>Description</b>         | groL antibody                                                        |
| <b>Species/Host</b>        | Rabbit                                                               |
| <b>Reactivity</b>          | Bacteria                                                             |
| <b>Conjugation</b>         | Unconjugated                                                         |
| <b>Tested Applications</b> | ELISA                                                                |
| <b>Immunogen</b>           | Recombinant Chlamydia pneumoniae 60 kDa chaperonin protein (1-214AA) |
| <b>Target</b>              | groL                                                                 |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol          |
| <b>Form/Appearance</b>     | Liquid                                                               |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.        |
| <b>Note</b>                | For research use only                                                |
| <b>Isotype</b>             | IgG                                                                  |
| <b>Clonality</b>           | Polyclonal                                                           |
| <b>Uniprot ID</b>          | <b>P31681</b>                                                        |
| <b>Dilution Range</b>      | ELISA                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                      |