

Product Datasheet

groL antibody (HRP) (orb51343)

Description	groL antibody (HRP)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Chlamydia pneumoniae 60 kDa chaperonin protein (1-214AA)
Target	groL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31681
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.