

Product Datasheet

EIF4B (phospho-Ser540) antibody (orb5137)

Description	EIF4B (phospho-Ser540) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human EIF4B around the phosphorylation site of Ser504 (QQ(p-S)PT)
Target	phospho-EIF4B (Ser504)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	69 kDa
Uniprot ID	P23588
Expiration Date	12 months from date of receipt.