

Product Datasheet

ompT antibody (orb51378)

Description	ompT antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli O157:H7 Protease 7 protein (21-317AA)
Target	ompT
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P58603
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.