

Product Datasheet

C11orf79 antibody (orb513800)

Description	Rabbit polyclonal antibody to C11orf79
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 1-166aa of human SDHAF2 protein
Target	SDHAF2
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	30KD
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.