

Product Datasheet

ompK antibody (FITC) (orb51384)

Description	ompK antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Vibrio parahaemolyticus serotype O3:K6 Outer membrane protein OmpK protein (21-266AA)
Target	ompK
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P59570
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.