

Product Datasheet

rpL antibody (FITC) (orb51396)

Description	rpL antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Staphylococcus aureus 50S ribosomal protein L7/L12 protein (1-122AA)
Target	rpL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P66061
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.