

Product Datasheet

IdhA antibody (Biotin) (orb51401)

Description	IdhA antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Staphylococcus aureus L-lactate dehydrogenase 1 protein (1-317AA)
Target	Idh1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P65255
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.