

Product Datasheet

Beta-lytic metalloendopeptidase antibody (Biotin) (orb51405)

Description	Beta-lytic metalloendopeptidase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Lysobacter enzymogenes Beta-lytic metalloendopeptidase protein (1-178AA)
Target	Beta-lytic metalloendopeptidase
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00801
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.