

Product Datasheet

eta antibody (FITC) (orb51416)

Description	eta antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Staphylococcus aureus Exfoliative toxin A protein (39-280AA)
Target	eta
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09331
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.