

Product Datasheet

gpd-1 antibody (Biotin) (orb51433)

Description	gpd-1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Other
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Caenorhabditis elegans Glyceraldehyde-3-phosphate dehydrogenase 1 protein (1-341AA)
Target	gpd-1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04970
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.