

Product Datasheet

lexA antibody (Biotin) (orb51449)

Description	lexA antibody (Biotin)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli LexA repressor protein (1-202AA)
Target	lexA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0A7C2
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.