

Product Datasheet

pyrE antibody (FITC) (orb51460)

Description	pyrE antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Laribacter hongkongensis Orotate phosphoribosyltransferase protein (1-213AA)
Target	pyrE
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	C1D6F5
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.