

Product Datasheet

IplA antibody (FITC) (orb51464)

Description	IplA antibody (FITC)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	FITC
Immunogen	Recombinant Escherichia coli Lipoate-protein ligase A protein (1-338AA)
Target	IplA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	C4ZT68
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.