


Product Datasheet

AK antibody (orb51466)

Description	AK antibody
Species/Host	Rabbit
Reactivity	Crustacean
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Penaeus monodon Arginine kinase protein (2-356AA)
Target	AK
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	C7E3T4
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant pro...