

Product Datasheet

pmp20 antibody (FITC) (orb51480)

Description	pmp20 antibody (FITC)
Species/Host	Rabbit
Reactivity	Fungi
Conjugation	FITC
Immunogen	Recombinant Neosartorya fumigata Peroxiredoxin Asp f3 protein (1-168AA)
Target	aspf3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43099
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.