

Product Datasheet

CLV1 receptor kinase-like protein antibody [Out of stock] (orb515153)

Description	Rabbit polyclonal antibody to CLV1 receptor kinase-like protein
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	A peptide corresponding to the C terminal of oryza sativa CLV1 receptor kinase-like protein.
Form/Appearance	Lyophilized power with 0.02% NaN ₃
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.