

Product Datasheet

Aquaporin (Water channel 1) (BfWC1) antibody [Out of stock] (orb515171)

Description	Rabbit polyclonal antibody to Aquaporin (Water channel 1) (BfWC1)
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	A peptide corresponding to the C terminal of haematobia irritans exigua Aquaporin protein.
Form/Appearance	Lyophilized power with 1% rAlbumin and 0.02% NaN ₃ .
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.