



Product Datasheet

ATG4D antibody (orb515194)

Description	Rabbit polyclonal antibody to ATG4D
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human ATG4D. The exact sequence is proprietary.
Target	ATG4D
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8BGV9 , Q86TLO
Entrez	84971 , 235040
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-100, IP: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of ATG4D expression...



Immunohistochemical analysis of ATG4D staining...