


Product Datasheet

TGM3 antibody (orb515217)

Description	Rabbit polyclonal antibody to TGM3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human Transglutaminase 3. The exact sequence is proprietary.
Target	TGM3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q08189 , D4A5U3 , Q08188
Entrez	7053 , 21818 , 366189
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Transglutaminase 3 (TGM3) in four lanes (A, B, C, D). Molecular weight markers are indicated on the left at 95, 72, 55, and 43 kDa.