



Product Datasheet

MKNK2 (phospho-T249) antibody (orb515263)


Description	Rabbit polyclonal antibody to MKNK2 (phospho-T249)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T249 of human MNK2 protein. The exact sequence is proprietary.
Target	MKNK2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q5U2N4 , Q8CDB0 , Q9HBH9
Entrez	17347 , 2872 , 299618
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of MNK2 (Phospho-T249)



Immunohistochemical analysis of MNK2 (Phospho-T249)



Direct ELISA antibody dose-response curve