



Product Datasheet

NOS2 (phospho-Y151) antibody (orb515265)


Description	Rabbit polyclonal antibody to NOS2 (phospho-Y151)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y151 of human iNOS protein. The exact sequence is proprietary.
Target	NOS2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P35228, Q06518, P29477
Entrez	4843, 18126, 24599
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of iNOS (Phospho-Y151)



Immunohistochemical analysis of iNOS (Phospho-Y151)



Direct ELISA antibody dose-response curve