



Product Datasheet

OTUB1 (phospho-S187) antibody (orb515276)


Description	Rabbit polyclonal antibody to OTUB1 (phospho-S187)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S187 of human OTUB1 protein. The exact sequence is proprietary.
Target	OTUB1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	B2RYG6, Q7TQI3, Q96FW1
Entrez	293705, 107260, 55611
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of OTUB1 (Phospho-S187) antibody.



Immunohistochemical analysis of OTUB1 (Phospho-S187) antibody.



Direct ELISA antibody dose-response curve for OTUB1 (Phospho-S187) antibody.