



Product Datasheet

ACIN1 (phospho-S1180) antibody (orb515284)

Description	Rabbit polyclonal antibody to ACIN1 (phospho-S1180)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S1180 of human Acinus protein. The exact sequence is proprietary.
Target	ACIN1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q14498, Q63627, O95104, Q8VH51, A2AJT4, Q8TF01, Q9JIX8, Q9UKV3
Entrez	57466, 56215, 22985, 25957, 9584, 170791
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Acinus (Phospho-S1180) antibody.



Direct ELISA antibody dose-response curve for Phosphopeptide and Non-phosphopeptide.