



Product Datasheet

PTK6 antibody (orb51542)

Description	PTK6 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Protein-tyrosine kinase 6 protein (1-113AA)
Target	PTK6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13882
Dilution Range	WB:1:2000-10000, IF:1:50-500
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 (lane 1), ...



Immunofluorescence analysis of MCF-7 cel...