



Product Datasheet

COL4A1 antibody (orb515470)

Description	Rabbit polyclonal antibody to COL4A1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human Collagen 4. The exact sequence is proprietary.
Target	COL4A1; COL4A2; COL4A4; COL4A5; COL4A6
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P53420 , Q14031 , P29400 , P08572 , P02462 , P02463
Entrez	12826 , 1282
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Collagen 4 expr...



Immunohistochemical analysis of Collagen...