



Product Datasheet

CEBPA (phospho-S21) antibody (orb515486)

Description	Rabbit polyclonal antibody to CEBPA (phospho-S21)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S21 of human C/EBP alpha protein. The exact sequence is proprietary.
Target	CEBPA
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P49715, P53566, P05554
Entrez	1050, 12606, 24252
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of C/EBP alpha (Phospho-S21)



Immunohistochemical analysis of C/EBP alpha (Phospho-S21)