



Product Datasheet

MITF (phospho-S180) antibody (orb515489)


Description	Rabbit polyclonal antibody to MITF (phospho-S180)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S180 of human MITF protein. The exact sequence is proprietary.
Target	MITF
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q08874, O75030
Entrez	17342, 4286
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of MITF (Phospho-S180) antibody.



Immunohistochemical analysis of MITF (Phospho-S180) antibody.



Direct ELISA antibody dose-response curve for MITF (Phospho-S180) antibody.