



Product Datasheet

CEBPE (phospho-T74) antibody (orb515496)

Description	Rabbit polyclonal antibody to CEBPE (phospho-T74)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T74 of human C/EBP epsilon protein. The exact sequence is proprietary.
Target	CEBPE
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q6PZD9 , P56261 , Q15744
Entrez	25410 , 110794 , 1053
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-100, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of C/EBP epsilon (...)



Immunohistochemical analysis of C/EBP epsilon...



Immunofluorescent analysis of C/EBP epsilon...