



Product Datasheet

FOSB (phospho-S27) antibody (orb515499)


Description	Rabbit polyclonal antibody to FOSB (phospho-S27)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S27 of human FOSB protein. The exact sequence is proprietary.
Target	FOSB
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P53539, P13346
Entrez	2354, 14282
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-100, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of FOSB (Phospho-S27) antibody.



Immunohistochemical analysis of FOSB (Phospho-S27) antibody.



Immunofluorescent analysis of FOSB (Phospho-S27) antibody.