



Product Datasheet

IGFBP3 (phospho-S183) antibody (orb515503)

Description	Rabbit polyclonal antibody to IGFBP3 (phospho-S183)
Species/Host	Rabbit
Reactivity	Human, Porcine
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S183 of human IGFBP3 protein. The exact sequence is proprietary.
Target	IGFBP3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P17936
Entrez	3486
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of IGFBP3 (Phospho-S183) antibody



Immunohistochemical analysis of IGFBP3 (Phospho-S183) antibody