



Product Datasheet

UBA52 (phospho-S65) antibody (orb515510)

Description	Rabbit polyclonal antibody to UBA52 (phospho-S65)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S65 of human Ubiquitin protein. The exact sequence is proprietary.
Target	UBA52; RPS27A; UBB; UBC
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P62987, P0CG47, P62979, P0CG48
Entrez	7311, 7314, 7316, 6233
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Ubiquitin (Phos...



Direct ELISA antibody dose-response curve...