



Product Datasheet

PRKCG (phospho-T514) antibody (orb515529)


Description	Rabbit polyclonal antibody to PRKCG (phospho-T514)
Species/Host	Rabbit
Reactivity	Human, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T514 of human PKC gamma protein. The exact sequence is proprietary.
Target	PRKCG
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P63319, P05129
Entrez	5582, 24681
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of PKC gamma (Phos...



Immunohistochemical analysis of PKC gamm...



Direct ELISA antibody dose-response curv...