



Product Datasheet

CDH1 (phospho-S844) antibody (orb515532)

Description	Rabbit polyclonal antibody to CDH1 (phospho-S844)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S844 of human E Cadherin protein. The exact sequence is proprietary.
Target	CDH1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9ROT4, P12830, P09803
Entrez	12550, 999, 83502
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of E Cadherin (Phospho-S844) antibody.



Direct ELISA antibody dose-response curve for Phosphopeptide and Non-phosphopeptide.