



Product Datasheet

PRKAG1 antibody (orb515562)

Description	Rabbit polyclonal antibody to PRKAG1
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human AMPK gamma 1/2/3. The exact sequence is proprietary.
Target	PRKAG1; PRKAG2; PRKAG3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9UGI9 , P54619 , O54950 , Q91WG5 , Q9UGJ0 , Q8BGM7 , P80385
Entrez	108099 , 25520 , 5571 , 53632 , 241113 , 51422 , 19082
Dilution Range	WB: 1:500-1000, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of AMPK gamma 1/2/...



Immunofluorescent analysis of AMPK gamma...