


Product Datasheet

GSK3B antibody (orb515586)

Description	Rabbit polyclonal antibody to GSK3B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant full length protein of human GSK3 beta. The exact sequence is proprietary.
Target	GSK3B
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P49841, P18266, Q9WV60
Entrez	84027, 2932, 56637
Dilution Range	WB: 1:500-1000, IP: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of GSK3 beta expression...