

Product Datasheet

GPR98 antibody (FITC) (orb515880)

Description	GPR98 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human G-protein coupled receptor 98 protein (6155-6306AA)
Target	GPR98
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WXG9
Expiration Date	12 months from date of receipt.