

Product Datasheet

ANO6 antibody (HRP) (orb515933)

Description	ANO6 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Anoctamin-6 protein (12-100AA)
Target	ANO6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q4KMQ2
Expiration Date	12 months from date of receipt.