



## Product Datasheet

### SCARB2 antibody (orb51610)


|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | SCARB2 antibody                                                  |
| <b>Species/Host</b>        | Rabbit                                                           |
| <b>Reactivity</b>          | Human, Mouse                                                     |
| <b>Conjugation</b>         | Unconjugated                                                     |
| <b>Tested Applications</b> | ELISA, IF, IHC, WB                                               |
| <b>Immunogen</b>           | Recombinant Human Lysosome membrane protein 2 protein (47-128AA) |
| <b>Target</b>              | SCARB2                                                           |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol      |
| <b>Form/Appearance</b>     | Liquid                                                           |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.    |
| <b>Note</b>                | For research use only                                            |
| <b>Isotype</b>             | IgG                                                              |
| <b>Clonality</b>           | Polyclonal                                                       |
| <b>Uniprot ID</b>          | <a href="#">Q14108</a>                                           |
| <b>Dilution Range</b>      | WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                  |



Western blot analysis of HepG2(lane 1), ...



Immunohistochemical staining of human co...



Immunofluorescence analysis of HepG2 cel...