

Product Datasheet

CNNM2 antibody (FITC) (orb516168)

Description	CNNM2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Metal transporter CNNM2 protein (739-875AA)
Target	CNNM2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H8M5
Expiration Date	12 months from date of receipt.