

Product Datasheet

SUN2 antibody (Biotin) (orb516214)

Description	SUN2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human SUN domain-containing protein 2 protein (45-155AA)
Target	SUN2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UH99
Expiration Date	12 months from date of receipt.